

First Part: Test Design

- 1- Verify that the user added an account from the accounts section
- 2- Verify that the user can delete an added account
- 3- Verify that the user can add a new transfer
- 4- Verify that the user can cancel the added transfer
- 5-Verify that the user adds new income
- 6-Verify that the user can add a new source of income
- 7-Verify that the user can change the currency of the app
- 8-Verify that the user can edit the category
- 9-Verify that the user can delete the category
- 10-Verify that the user added a new expense

Second Part: Bug Reporting

1-

Title → password filed automatically when user sign up with mail

Reproducible Steps →

user installs the app

Press login button

Insert the phone number


Insert the received OTP

Fulfill the field

Insert the email

Don't insert data in the password field

Attachments →



Affected Devices → oppo reno 2f device Android 11


Network → 5G

Severity → medium
Priority → medium
Impact → the password filed Fulfill automatically
Type → Functionality bug

2-

Title → All notifications and blog sections appear in Arabic while the app language is English
Reproducible Steps →
Choose the English language for the app
Press on the side menu
Select Blog
Check the language

Attachments →





Affected Devices → oppo reno 2f device Android 11
Network → 5G
Severity → medium
Priority → high
Impact → the language in the blog appears in Arabic and all notifications of the app are sent in Arabic
Type → localization bug

3-

Title → The wishlist is empty when the user adds a favorite product
Reproducible Steps →
User login
Add product to wishlist
Press on the side menu

Check the wishlist

Attachments →



Affected Devices → oppo reno 2f device Android 11

Network → 5G
Severity → high
Priority → high
Impact → The product must be added to the wishlist, but the wishlist is empty.
type → Functionality bug

4-

Title → There is an incorrect number that appears when there is an increase or decrease in the product quantity

Reproducible Steps →

User login

The user increases the quantity of any product

The user decreases the quantity of the same product

And make it many times

Attachments →

The bug appears in the screen record

<https://drive.google.com/file/d/18OOt5pdoESB5hXCfaBfjTwGUGo3lrjF/view?usp=sharing>

Affected Devices → oppo reno 2f device Android 11

Network → 5G

Severity → critical

Priority → high

Impact → incorrect number that appears when there is an increase or decrease in the product quantity.

Type → Functionality Bug

5-

Title → Numbers of height and weight were written incorrectly when editing the health profile.

Reproducible Steps →

User login

Press on the side menu

Select health profile

Insert data in height and weight

Check the inserted data

Attachments →

The screenshot shows a mobile app interface for editing a 'Health profile'. At the top, the status bar shows the time 3:06 and various icons. Below the title 'Health profile', there are two input fields: 'Weight' with the value '0.060' and 'Height' with the value '0.0160'. Both fields are circled in red. Below these is a 'Blood type' dropdown menu set to 'O+'. Under the 'Chronic Diseases' section, there is a green bar indicating 'Hypertension (High blood pressure)' is checked. A green 'Save' button is located below the diseases section. At the bottom, there is a numeric keypad with digits 1-9, 0, a comma, and a decimal point, along with navigation icons.

Affected Devices → oppo reno 2f device Android 11

Network → 5G

Severity → medium

Priority → medium

Impact → Numbers of height and weight were written incorrectly when editing the health profile.

Type → usability bug